Synovis®

Surgical Innovations

INSTRUCTIONS FOR USE	2
MODE D'EMPLOI	5
GEBRAUCHSANLEITUNG	8
ISTRUZIONI PER L'USO	11
INSTRUCCIONES DE USO	14
GEbruiksinstructies	17
BETJENINGSVEJLEDNING	20
ANVÄNDARINSTRUKTIONER	23
BRUKSANVISNINGER	26
KULLANMA TALIMATI	29
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	32
KÄYTTÖOHJEET	36
INSTRUÇÕES DE UTILIZAÇÃO	39
INSTRUCȚIUNI DE FOLOSIRE	42
NÁVOD K POUŽITÍ	45
HASZNÁLATI ÚTMUTATÓ	48
INSTRUKCJA UŻYCIA	51

Dura - Guard®

Dural Repair Patch • with Apex Processing®
Patch pour réparation durale • avec Apex Processing®
Dura-Reparaturpatch • mit Apex Processing®
Patch per riparazioni della dura madre • con Apex Processing®
Parche de reparación durale • con Apex Processing®
Pleisters voor herniaherstel • met Apex Processing®
Dural rekonstruktionsstykker • med Apex Processing®
Dural reparationspatch • med Apex Processing®
Dura mater rekonstruksjonslapp • med Apex Processing®
Dural Onarım Yaması • ile Apex Processing®
Εμβάλωμα αποκατάστασης σκληράς μήνιγγας • με Apex Processing®
Paikkamateriaali duuran korjauksiin • mukana Apex Processing®
Retalho para Reparação da dura-máter • com Apex Processing®
Bandă pentru repararea defectelor durale • cu Apex Processing®
Záplata pro reparaci dury • s technologií zpracování Apex Processing®
Dura-helyreállító tapasz • Apex Processing®-gel
Łata do rekonstrukcji opony twardej • z Apex Processing®

CE
2797

SYMBOL DEFINITIONS:

DO NOT FREEZE

Do not freeze

0°C / 32°F

Lower limit of temperature

Do not reuse

Consult *Instructions for Use*

STERILE A

Sterilized using aseptic processing techniques

NaOH

This product is treated with sodium hydroxide

BOVINE

This product is derived from U.S.D.A.-inspected cattle

MADE IN THE U.S.A.

Made in the U.S.A.

Rx Only

CAUTION: Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician

CONTENT

Content

REF

Catalog number

Use by date

LOT

Lot number

Manufacturer

EC REP

Authorized Representative in the European Community

PN

SSI part number

IC

SSI internal code

INDICATIONS:

For use as a dura substitute for the closure of dura mater during neurosurgery.

DESCRIPTION:

Dura-Guard® is prepared from bovine pericardium which is cross-linked with glutaraldehyde. The pericardium is procured from cattle originating in the United States. Dura-Guard is chemically sterilized using ethanol and propylene oxide. Dura-Guard has been treated with 1 molar sodium hydroxide for 60-75 minutes at 20-25°C.

Dura-Guard is packaged in a container filled with sterile, non-pyrogenic water containing propylene oxide. The contents of the unopened, undamaged container are sterile.

Dura-Guard utilizes animal tissue; patient must be informed prior to any procedure.

CONTRAINDICATIONS:

Dura-Guard is not designed, sold or intended for use except as indicated.

INSTRUCTIONS FOR USE:

I. RINSE PROCEDURE

1. Remove inner container from outer cardboard package. Do not place the container in the sterile field.

2. Examine the freeze indicator. Do not use if activated.
3. Inspect container and package. Do not use if there is evidence of moisture or leakage.
4. Open inner container. Use sterile, atraumatic forceps to grasp the edge of the patch and remove from the container using aseptic technique. Rinse surgical gloves to remove glove powder prior to touching the patch.
5. Immerse and agitate Dura-Guard, for a minimum of 3 minutes, in a sterile basin containing 500 ml of sterile physiologic saline (0.9% NaCl). Do not pour the storage solution into the sterile physiologic saline.

At the surgeon's discretion the 500 ml rinse solution may contain one of the following antibiotic treatments: ampicillin & gentamicin, bacitracin, cefazolin, cefotaxime, neomycin, and vancomycin. Testing has shown that Dura-Guard is not adversely affected by treatment with the antibiotics listed. The effects of other antibiotics on Dura-Guard have not been tested. The long term effects of antibiotic treatments of Dura-Guard have not been assessed. Do not use antibiotics contrary to the antibiotic manufacturer's instructions.

6. Keep the patch immersed in sterile saline until ready to use.
THE PATCH MUST REMAIN MOIST AT ALL TIMES.

II. IMPLANT INSTRUCTIONS

1. The patch may be tailored during surgery to meet the surgeon's needs.
2. Visually examine both sides of the Dura-Guard patch. If one side appears smoother, implant the smoother surface so that it faces the neural tissue.
3. Dura-Guard may be sutured, clipped, or stapled to the edge of the host tissue.

4. When implanting by suture, suture bites should be taken 2 to 3 millimeters from the edge of the graft.
5. The graft should be applied and fixed in place carefully to obtain best results.
6. Discard any unused pieces of Dura-Guard. Dura-Guard cannot be resterilized or reused.

STORAGE:

Do not freeze. Store at room temperature.

WARNINGS:

The Rinse Procedure must be followed or a sterile inflammatory reaction in the adjoining host tissue may result. Do not pour the storage solution into the rinse bath.

To avoid damage to the product, do not expose to any chemicals or substances other than those specified in this Rinse Procedure. Antimycotics must not come in contact with Dura-Guard as they are believed to alter the cross-link characteristics of tissue fixed in aldehyde preparations.

Do not freeze. Damage may result. Do not use if freeze indicator is activated.

Do not resterilize. Do not subject to steam, gas (ethylene oxide), or radiation sterilization as these may damage Dura-Guard.

Do not use if the container is not properly sealed prior to opening, as sterility may be compromised. Do not place container in sterile field as the outside of the container is not sterile. This product is for single use only. Once the container seal is broken, use immediately, do not reseal container or reuse Dura-Guard. Any unused pieces of Dura-Guard must be discarded. Failure to observe these warnings may result in surgical infection.

PRECAUTIONS:

Rinse surgical gloves to remove glove powder prior to touching Dura-Guard.

Before surgery, prospective patients or their representatives should be informed about complications which may be associated with use of this product. As with any surgical procedure, dehiscence at the surgical site and infection are possible complications. Monitor patient for infection and take appropriate therapeutic action.

ADVERSE REACTIONS:

Failure to rinse the product may result in a sterile inflammatory reaction (see Warnings and Rinse Procedure).

The effect on long term surgical outcome of damaging bovine pericardium by contact with chemicals or substances (other than saline), by freezing, or by exposure to steam, gas (ethylene oxide), or radiation sterilization has not been investigated (see Warnings).

Use of product following compromise of sterility may result in infection (see Warnings).

Long term incidence rates of host reactions (calcification, infection, rejection, adhesion, and hematological compatibility) during use for dural repair have not been investigated.

DISCLAIMER OF WARRANTIES:

Synovis Surgical Innovations (SSI), a division of Synovis Life Technologies, Inc., warrants that reasonable care has been used in the manufacture of this device. This warranty is exclusive and in lieu of all other warranties whether expressed, implied, written or oral, including, but not limited to, any implied warranties of merchantability or fitness. As a result of biological differences in individuals, no product is 100% effective under all circumstances. Because of this fact and since SSI has no control over the conditions under which the device is used, diagnosis of the patient, methods of administration or its handling after it leaves its possession, SSI does not warrant either a good effect or against an ill effect following its use. The manufacturer shall not be liable for any incidental or consequential loss, damage or expense arising directly or indirectly from the use of this device. SSI will replace any device which is defective at the time of shipment. No representative of SSI may change any of the foregoing or assume any additional liability or responsibility in connection with this device.

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Dura-Guard® is a registered trademark of Synovis Life Technologies, Inc.
A Division of Synovis Life Technologies, Inc.