
<b>(ENG)</b>	Instructions for Use	<b>(DUT)</b>	Gebruiksaanwijzing
<b>(CZE)</b>	Návod k použití	<b>(NOR)</b>	Bruksanvisning
<b>(DAN)</b>	Brugsanvisning	<b>(POL)</b>	Instrukcja Stosowania
<b>(GER)</b>	Gebrauchsanweisung	<b>(POR)</b>	Instruções de Utilização
<b>(GRE)</b>	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	<b>(SWE)</b>	Instruktioner för Användning
<b>(SPA)</b>	Instrucciones de Uso	<b>(BUL)</b>	Инструкция за употреба
<b>(FIN)</b>	Käyttöohje	<b>(ROU)</b>	Instrucțiune de utilizare
<b>(FRE)</b>	Instructions d'utilisation	<b>(HRV)</b>	Upute za uporabu
<b>(ITA)</b>	Istruzioni per l'uso	<b>(TUR)</b>	Kullanma Talimatları
<b>(LIT)</b>	Naudojimo instrukcija	<b>(BRA)</b>	Instruções de Utilização

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
<p><b>DESCRIPTION</b></p> <p>OSTENE is a water-soluble surgical implant material, and it can be used for the control of bleeding from bone surfaces by functioning as a mechanical barrier. OSTENE is a sterile mixture of water-soluble alkylene oxide copolymers.</p>
<p><b>INDICATIONS</b></p> <p>OSTENE is indicated for use as a water-soluble implant material and for use in the control of bleeding from bone surfaces. OSTENE should be warmed to desired consistency using aseptic techniques.</p>
<p><b>ACTIONS</b></p> <p>OSTENE achieves local hemostasis of bone by functioning as a mechanical barrier. It does not act biochemically.</p>
<p><b>CONTRAINDICATIONS</b></p> <ol style="list-style-type: none"><li>Sites with active or latent infections.</li> <li>The use of the product to lend structural support to bone.</li></ol>
<p><b>WARNINGS AND PRECAUTIONS</b></p> <ol style="list-style-type: none"><li>OSTENE is provided sterile for single use only. DO NOT RESTERILIZE. Discard any open unused, or damaged packages. DO NOT USE if there is a loss of sterility of the device. Once the sterile boundary is broken the device cannot be re-used because there is a risk of contamination.</li> <li>OSTENE should not be subjected to excess heat. Store at 0°C to 32°C (32°F to 90°F), away from direct heat, including sunlight.</li></ol>

**INSTRUCTIONS FOR USE**

OSTENE works best at body temperature and should be warmed to desired consistency using aseptic technique and then manipulated and softened with dry gloved fingers. The unopened foil packet may be immersed into a warm sterile solution at temperatures at or below 43°C (110°F). OSTENE handles more like a gum or putty than a wax. OSTENE may be pressed into bleeding bone and may be worked across the bleeding surface as indicated by the surgical circumstances and the preference of the surgeon. When using in the control of bone bleeding following a sternotomy, one packet or more may be applied to each edge of the cut sternum.

**STERILITY**

OSTENE is sterilized using electron beam irradiation. Removal from the sterile package using aseptic techniques should only take place immediately before use. OSTENE MUST NOT BE RESTERILIZED BY ANY METHOD.

**WARRANTY**

All products are warranted to be free from defects in material and workmanship. No warranty is made for any purpose other than in the product specifications and labeling.

**CAUTION**

Federal law (USA) restricts this device to sale by or on the order of a physician.

<p><b>(ENG)</b></p> <p>Instructions for Use</p>	<p><b>(DUT)</b></p> <p>Gebruiksaanwijzing</p>
<p><b>(CZE)</b></p> <p>Návod k použití</p>	<p><b>(NOR)</b></p> <p>Bruksanvisning</p>
<p><b>(DAN)</b></p> <p>Brugsanvisning</p>	<p><b>(POL)</b></p> <p>Instrukcja Stosowania</p>
<p><b>(GER)</b></p> <p>Gebrauchsanweisung</p>	<p><b>(POR)</b></p> <p>Instruções de Utilização</p>
<p><b>(GRE)</b></p> <p>ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ</p>	<p><b>(SWE)</b></p> <p>Instruktioner för Användning</p>
<p><b>(SPA)</b></p> <p>Instrucciones de Uso</p>	<p><b>(BUL)</b></p> <p>Инструкция за употреба</p>
<p><b>(FIN)</b></p> <p>Käyttöohje</p>	<p><b>(ROU)</b></p> <p>Instrucțiune de utilizare</p>
<p><b>(FRE)</b></p> <p>Instructions d'utilisation</p>	<p><b>(HRV)</b></p> <p>Upute za uporabu</p>
<p><b>(ITA)</b></p> <p>Istruzioni per l'uso</p>	<p><b>(TUR)</b></p> <p>Kullanma Talimatları</p>
<p><b>(LIT)</b></p> <p>Naudojimo instrukcija</p>	<p><b>(BRA)</b></p> <p>Instruções de Utilização</p>

**OSTENE** is a water-soluble surgical implant material, and it can be used for the control of bleeding from bone surfaces by functioning as a mechanical barrier. OSTENE is a sterile mixture of water-soluble alkylene oxide copolymers.

## OSTENE

**DESCRIPTION**

OSTENE is a water-soluble surgical implant material, and it can be used for the control of bleeding from bone surfaces by functioning as a mechanical barrier. OSTENE is a sterile mixture of water-soluble alkylene oxide copolymers.

**INDICATIONS**

OSTENE is indicated for use as a water-soluble implant material and for use in the control of bleeding from bone surfaces. OSTENE should be warmed to desired consistency using aseptic techniques.

**ACTIONS**

OSTENE achieves local hemostasis of bone by functioning as a mechanical barrier. It does not act biochemically.

**CONTRAINDICATIONS**

- Sites with active or latent infections.
- The use of the product to lend structural support to bone.

<p><b>(ENG)</b></p> <p>Instructions for Use</p>	<p><b>(DUT)</b></p> <p>Gebruiksaanwijzing</p>
<p><b>(CZE)</b></p> <p>Návod k použití</p>	<p><b>(NOR)</b></p> <p>Bruksanvisning</p>
<p><b>(DAN)</b></p> <p>Brugsanvisning</p>	<p><b>(POL)</b></p> <p>Instrukcja Stosowania</p>
<p><b>(GER)</b></p> <p>Gebrauchsanweisung</p>	<p><b>(POR)</b></p> <p>Instruções de Utilização</p>
<p><b>(GRE)</b></p> <p>ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ</p>	<p><b>(SWE)</b></p> <p>Instruktioner för Användning</p>
<p><b>(SPA)</b></p> <p>Instrucciones de Uso</p>	<p><b>(BUL)</b></p> <p>Инструкция за употреба</p>
<p><b>(FIN)</b></p> <p>Käyttöohje</p>	<p><b>(ROU)</b></p> <p>Instrucțiune de utilizare</p>
<p><b>(FRE)</b></p> <p>Instructions d'utilisation</p>	<p><b>(HRV)</b></p> <p>Upute za uporabu</p>
<p><b>(ITA)</b></p> <p>Istruzioni per l'uso</p>	<p><b>(TUR)</b></p> <p>Kullanma Talimatları</p>
<p><b>(LIT)</b></p> <p>Naudojimo instrukcija</p>	<p><b>(BRA)</b></p> <p>Instruções de Utilização</p>

<p><b>(ENG)</b></p> <p>Instructions for Use</p>	<p><b>(DUT)</b></p> <p>Gebruiksaanwijzing</p>
<p><b>(CZE)</b></p> <p>Návod k použití</p>	<p><b>(NOR)</b></p> <p>Bruksanvisning</p>
<p><b>(DAN)</b></p> <p>Brugsanvisning</p>	<p><b>(POL)</b></p> <p>Instrukcja Stosowania</p>
<p><b>(GER)</b></p> <p>Gebrauchsanweisung</p>	<p><b>(POR)</b></p> <p>Instruções de Utilização</p>
<p><b>(GRE)</b></p> <p>ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ</p>	<p><b>(SWE)</b></p> <p>Instruktioner för Användning</p>
<p><b>(SPA)</b></p> <p>Instrucciones de Uso</p>	<p><b>(BUL)</b></p> <p>Инструкция за употреба</p>
<p><b>(FIN)</b></p> <p>Käyttöohje</p>	<p><b>(ROU)</b></p> <p>Instrucțiune de utilizare</p>
<p><b>(FRE)</b></p> <p>Instructions d'utilisation</p>	<p><b>(HRV)</b></p> <p>Upute za uporabu</p>
<p><b>(ITA)</b></p> <p>Istruzioni per l'uso</p>	<p><b>(TUR)</b></p> <p>Kullanma Talimatları</p>
<p><b>(LIT)</b></p> <p>Naudojimo instrukcija</p>	<p><b>(BRA)</b></p> <p>Instruções de Utilização</p>

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