Synovis®

Surgical Innovations

INSTRUCTIONS FOR USE	2
MODE D'EMPLOI.....	4
ANWENDUNGSHINWEISE	6
ISTRUZIONI PER L'USO.....	8
INSTRUCCIONES DE USO.....	10
IGEBRUIKSAANWIJZING.....	12
IBRUGSANVISNING	14
BRUKSANVISNING	16
BRUKSANVISNING	18
KULLANMA TALIMATI.....	20
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ.....	22
KÄYTTÖOHJEET	24
INSTRUÇÕES DE UTILIZAÇÃO.....	26
ISTRUCȚIUNI DE FOLOSIRE.....	28
POKYNY K POUŽITÍ.....	30
HASZNÁLATI ÚTMUTATÓ	32
INSTRUKCJA STOSOWANIA	34

PSD GEL

Adhesive for Peri-Strips Dry®
Adhésif pour Peri-Strips Dry®
Klebstoff für Peri-Strips Dry®
Adesivo per le strisce Peri-Strips Dry®
Adhesivo para Peri-Strips Dry®
Hechtmateriaal voor Peri-Strips Dry®
Bindemiddel til Peri-Strips Dry®
Adhesiv för Peri-Strips Dry®
Lim for Peri-Strips Dry®
Peri-Strips Dry® için Yapışkan
Αυτοκόλλητο για Peri-Strips Dry®
Peri-Strips Dry®:lle tarkoitettu sidosaine
Adesivo para Peri-Strips Dry®
Adeziv pentru Peri-Strips Dry®
Lepidlo pro Peri-Strips Dry®
Ragasztóanyag a Peri-Strips Dry® termékhez
Klej do Peri-Strips Dry®

SYMBOL DEFINITIONS:

Do not reuse

Consult *Instructions for Use*

STERILE

Sterile

STERILE **EO**

Gel package is Ethylene Oxide sterilized.

STERILE

Gel is steam sterilized.

MADE IN THE U.S.A.

Made in the U.S.A.

Rx Only

CAUTION: Federal (USA) law restricts this device to sale by, or on the order of, a physician.

CONTENT

Content

REF

Catalog number

Use by date

LOT

Lot number

Manufacturer

EC **REP**

Authorized Representative in the European Community

PN

SSI part number

TN

SSI tracking number

IC

SSI internal code

DESCRIPTION:

PSD Gel is prepared from Hydroxypropyl Methyl Cellulose, Propylene Glycol, and water.

The PSD Gel is used to bond the Peri-Strips Dry pericardium strips onto the surgical stapler forks until the stapler is positioned and fired. Each PSD Gel tube is packaged sterile in a pouch.

INDICATIONS:

PSD Gel is intended to be used only with Peri-Strips Dry® (PSD) as a temporary bonding agent during applications indicated in the Peri-Strips Dry Instructions for Use.

CONTRAINDICATIONS:

PSD Gel is not designed, sold or intended for use except as indicated. Use only as a temporary bonding agent with Peri-Strips Dry.

WARNINGS/PRECAUTIONS:

The contents of the unopened, undamaged pouches are sterile. Do not re-sterilize PSD Gel.

Do not use if pouches are not properly sealed prior to opening, as sterility may be compromised. Do not place pouches in the sterile field, as the outside of the pouches are not sterile. This product is for single use only. Once any pouch seal is broken, use immediately. Any opened pouches of PSD Gel must be discarded.

Failure to observe these warnings may result in surgical infection.

Do not get the Peri-Strips Dry® pericardium strips wet prior to applying PSD Gel, as this may prevent bonding. Keep the surgical site as dry as possible to prevent the Peri-Strips Dry® pericardium strips from shifting or falling off the stapler.

INSTRUCTIONS FOR USE:

1. Inspect the pouch. Do not use if there is evidence of punctures or tears.
2. Peel open the pouch and aseptically remove the PSD Gel tube; the PSD Gel tube may be placed in the sterile field.
3. Break off the metal tip of the PSD Gel tube.
4. Attach the plastic gel tip to the top of the PSD Gel tube.
5. Ensure the anvil and cartridge side of the stapler is clean and dry.
6. Apply a thin continuous line of PSD Gel following the instructions given in the PSD staple line reinforcement product IFU for your specific device.
7. Discard any opened PSD Gel tubes.

STORAGE:

Store PSD Gel at room temperature (15°-30°C / 59°-86°F).

DISCLAIMER OF WARRANTIES:

Synovis Surgical Innovations (SSI), a division of Synovis Life Technologies, Inc., warrants that reasonable care has been used in the manufacture of this device. This warranty is exclusive and in lieu of all other warranties whether expressed, implied, written or oral, including, but not limited to, any implied warranties of merchantability or fitness. As a result of biological differences in individuals, no product is 100% effective under all circumstances. Because of this fact and since SSI has no control over the conditions under which the device is used, diagnosis of the patient, methods of administration or its handling after it leaves its possession, SSI does not warrant either a good effect or against an ill effect following its use. The manufacturer shall not be liable for any incidental or consequential loss, damage or expense arising directly or

indirectly from the use of this device. No representative of SSI may change any of the foregoing or assume any liability or responsibility in connection with this device.

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A Division of Synovis Life Technologies, Inc.
AR-MED, a division of Quintiles Ltd,