Baxter

DURA-GUARD

- (en) DURAL REPAIR PATCH
- (fr) PATCH DE RÉPARATION DE LA DURE-MÈRE
- (de) DURA-REPARATURPATCH
- (it) PATCH DI RIPARAZIONE DELLA DURA
- (es) PARCHE DE REPARACIÓN DE LA DURAMADRE
- (nl) DURALE HERSTELPATCH
- (da) DURAREPARATIONSPLASTER
- (sv) PLÅSTER FÖR DURALREPARATION
- (no) DURA-REPARASJONSLAPP
- (tr) DURAL ONARIM YAMASI
- (el) ΕΠΙΘΕΜΑ ΑΠΟΚΑΤΑΣΤΑΣΗΣ ΣΚΛΗΡΑΣ ΜΗΝΙΓΓΑΣ
- (fi) KOVAKALVON KORJAUSPAIKKA
- (pt) PATCH DE REPARAÇÃO DA DURA-MÁTER
- (ro) MATERIAL PENTRU REPARAȚII DURALE
- (cs) NÁPLAST PRO DURÁLNÍ OPRAVU
- (hu) DURA-HELYREÁLLÍTÓ FOLT
- (pl) ŁATA DO OPONY TWARDEJ

INSTRUCTIONS FOR USE (IFU).....	2
MODE D'EMPLOI	6
GEBRAUCHSANWEISUNG.....	10
ISTRUZIONI PER L'USO	14
INSTRUCCIONES DE USO	18
GEBRUIKSAANWIJZING	22
BRUGSANVISNING (IFU).....	26
BRUKSANVISNING	30
BRUKSANVISNING	34
KULLANIM TALİMATI (IFU).....	38
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	42
KÄYTTÖOHJEET.....	46
INSTRUÇÕES DE UTILIZAÇÃO	50
INSTRUCȚIUNI DE UTILIZARE (IDU).....	54
NÁVOD K POUŽITÍ	58
HASZNÁLATI UTASÍTÁS	62
INSTRUKCJA UŻYCIA.....	66

Symbols Referenced on Labeling

Symbol Glossary per US FD&C Act:

Standard	Symbol	Symbol Title	Symbol Meaning	Symbol Number
ISO 15223-1*		Manufacturer	Manufacturer	5.1.1
ISO 15223-1		Authorized representative in the European Community	Authorized representative in the European Community	5.1.2
ISO 15223-1		Use-by date	Use-by date	5.1.4
ISO 15223-1		Batch code	Lot number	5.1.5
ISO 15223-1		Catalogue number	Catalog number	5.1.6
ISO 15223-1		Sterilized using aseptic processing techniques	Sterilized using aseptic processing techniques	5.2.2
ISO 15223-1		Do not re-sterilize	Do not re-sterilize	5.2.6
ISO 15223-1		Do not use if package is damaged	Do not use if the product sterile barrier or its packaging is compromised	5.2.8

Standard	Symbol	Symbol Title	Symbol Meaning	Symbol Number
ISO 15223-1		Temperature limit	Store at controlled room temperature	5.3.7
ISO 15223-1		Do not re-use	Do not re-use	5.4.2
ISO 15223-1		Consult instructions for use	Consult instructions for use	5.4.3
ISO 15223-1		Caution	Caution: Consult instruction for use for warning and precaution information	5.4.4
			Content	
			Caution: Federal (USA) law restricts this device to sale by or on the order of a physician	

*ISO 15223-1: 2016, Medical Devices-Symbols to be used with medical device labels, labeling and information to be supplied, Part 1: General requirements

Additional symbols on the product labeling that are not required by the US FD&C Act:

Symbol	Symbol Meaning
MADE IN THE U.S.A.	Made in the U.S.A.
PN	Manufacturer part number
DO NOT FREEZE	Do not freeze
BOVINE	This product is derived from cattle
MODEL	Model number
SEE IFU FOR SYMBOL DEFINITIONS	See IFU for symbol definitions

DESCRIPTION:

DURA-GUARD is prepared from bovine pericardium which is cross-linked with glutaraldehyde. The pericardium is procured from cattle originating in the United States. DURA-GUARD is chemically sterilized using ethanol and propylene oxide. DURA-GUARD has been treated with 1 molar sodium hydroxide for a minimum of 60 minutes at 20 - 25 °C.

DURA-GUARD is packaged in a container filled with sterile, non-pyrogenic water containing propylene oxide. The contents of the unopened, undamaged container are sterile.

DURA-GUARD is MR Safe.

DURA-GUARD utilizes animal tissue; patient must be informed prior to any procedure.

DURA-GUARD Product Models and Sizes

DURA-GUARD Model Number	Size (cm)
DG-0209SN	2X9
DG-0404SN	4X4
DG-0608SN	6X8
DG-0814SN	8X14
DG-1016SN	10X16

INDICATIONS:

For use as a dura substitute for the closure of dura mater during neurosurgery.

CONTRAINDICATIONS:

DURA-GUARD is not designed, sold or intended for use except as indicated.

Do not use DURA-GUARD in patients with a known sensitivity to bovine material.

ADVERSE REACTIONS:

As with any surgical procedure, dehiscence at the surgical site, hematoma, seroma, fever and infection are possible complications. Other

potential complications specific to dura-related procedures include, but are not limited to: death, seizure, cerebrospinal fluid leakage or fistula, meningitis, adhesion formation, pseudomeningocele or hydrocephalus. Monitor patient for adverse reactions and take appropriate therapeutic action.

Failure to rinse the product may result in a sterile inflammatory reaction (see Warnings and Rinse Procedure).

The effect on long term surgical outcome of damaging bovine pericardium by contact with chemicals or substances (other than saline), by freezing, or by exposure to steam, gas (ethylene oxide), or radiation sterilization has not been investigated (see Warnings).

Use of product following compromise of sterility may result in infection (see Warnings).

Long term incidence rates of host reactions (calcification, infection, rejection, and adhesion) during use for dural repair have not been investigated.

Glutaraldehyde-treated bovine pericardium may undergo accelerated calcific infiltration in patients with high calcium metabolic activity (e.g., children). This may not be a concern where the patch is exposed to systolic pressures.

WARNINGS:

The Rinse Procedure must be followed or a sterile inflammatory reaction in the adjoining host tissue may result. Do not pour the storage solution into the rinse bath.

To avoid damage to the product, do not expose to any chemicals or substances other than those specified in this Rinse Procedure. Antimycotics (anti-fungals) must not come in contact with DURA-GUARD as they are believed to alter the cross-link characteristics of tissue fixed in aldehyde preparations.

Do not freeze. Damage to the sterile barrier may result. Do not use if freeze indicator is activated.

The product must remain moist at all times.

Do not use DURA-GUARD in patients with a known sensitivity to bovine material.

Synovis products differ; substitution of one product for another product may be harmful to the patient.

Do not resterilize. Do not subject to steam, gas (ethylene oxide), or radiation sterilization as these may damage DURA-GUARD.

Do not use if the container is not properly sealed prior to opening, as sterility may be compromised. Do not place container in sterile field as the outside of the container is not sterile. This product is for single use only. Once the container seal is broken, use immediately, do not reseal container or reuse DURA-GUARD. Any unused pieces of DURA-GUARD must be discarded as biohazardous waste. Failure to observe these warnings may result in surgical infection.

PRECAUTIONS:

Rinse surgical gloves to remove glove powder prior to touching DURA-GUARD.

Before surgery, prospective patients or their representatives should be informed about complications which may be associated with use of this product.

INSTRUCTIONS FOR USE:

I. RINSE PROCEDURE

1. Remove inner container from outer cardboard package. Do not place the container in the sterile field.
2. Examine the freeze indicator located inside the lid of the outer cardboard package. Do not use if activated.
3. Inspect container and package. Verify shrink wrap and tamper-evident tape are intact. Do not use if there is evidence of tampering, moisture or leakage.
4. Open inner container. Use sterile, atraumatic forceps to grasp the edge of the patch and remove from the container using aseptic technique. Rinse surgical gloves to remove glove powder prior to touching the patch.

5. Immerse and agitate DURA-GUARD, for a minimum of 3 minutes, in a sterile basin containing 500 ml of sterile physiologic saline (0.9% NaCl). Do not pour the storage solution into the sterile physiologic saline.

At the surgeon's discretion the 500 ml rinse solution may contain one of the following antibiotic treatments: ampicillin & gentamicin, bacitracin, cefazolin, cefotaxime, neomycin, and vancomycin. Testing has shown that DURA-GUARD is not adversely affected by treatment with the antibiotics listed. The effects of other antibiotics on DURA-GUARD have not been tested. The long term effects of antibiotic treatments on DURA-GUARD have not been assessed. Do not use antibiotics contrary to the antibiotic manufacturer's instructions.

6. Keep the patch immersed in sterile saline until ready to use. THE PATCH MUST REMAIN MOIST AT ALL TIMES.

II. IMPLANT INSTRUCTIONS

1. The patch may be tailored during surgery to meet the surgeon's needs.
2. DURA-GUARD may be sutured to the edge of the host tissue.
3. When implanting by suture, suture bites should be taken 2 to 3 millimeters from the edge of the graft.
4. The graft should be applied and fixed in place carefully to obtain best results.
5. Discard any unused pieces of DURA-GUARD in biohazardous waste. DURA-GUARD cannot be resterilized or reused.

STORAGE:

Do not freeze. Store at room temperature 20 °C – 25 °C (68 °F – 77 °F).

DISPOSAL

Any packaging or components exposed to human tissue/fluids should be disposed of per hospital protocols. Any opened, unused components should be discarded due to compromised sterility.

DISCLAIMER OF WARRANTIES:

Synovis Life Technologies, Inc. (SLT) (A Subsidiary of Baxter International Inc.) warrants that reasonable care has been used in the manufacture of this device. As a result of biological differences in individuals, no product is 100% effective under all circumstances. Because of this fact and since SLT has no control over the conditions under which the device is used, diagnosis of the patient, methods of administration or its handling after it leaves its possession, SLT does not warrant either a good effect or against an ill effect following its use. SLT will replace any device, which is defective at the time of shipment. No representative of SLT may change any of the foregoing or assume any additional liability or responsibility in connection with this device.

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Baxter and Dura-Guard are registered trademarks of Baxter International Inc., or its subsidiaries.