

Tisseel

FIBRIN SEALANT

INTRODUCING GASLESS SPRAY TIP

For use with TISSEEL Fibrin Sealant

- The new applicator technology that is ready when you are



Saves Time

Designed to reduce OR setup time by eliminating the need for external gas sources¹



Dual functionality

Gasless Spray Tip delivers your choice of broad spray or precise drip¹



Open or Laparoscopic

Your choice of Open, 25 cm MIS and 40 cm MIS spray applicators



Quick Connect

Secure Snap-Lock attachment to PRIMA syringe¹



Baxter

GASLESS SPRAY TIP

The new applicator technology that is ready when you are

Open and MIS Configurations

The Gasless Spray Tip Applicator is available in three configurations: Open, 25cm MIS, and 40cm MIS. Each applicator pack contains two replacement spray tips.



Gasless Spray Tip Applicator for Open procedures

ORDER CODE 601151 QTY 5 per CA



Gasless Spray Tip Applicator for MIS procedures, 25cm

ORDER CODE 601153 QTY 5 per CA



Gasless Spray Tip Applicator for MIS procedures, 40cm

ORDER CODE 601155 QTY 5 per CA

GASLESS TIP INDICATIONS FOR USE

The Gasless Tip is indicated for the application of two non-homogenous liquids.

WARNING/PRECAUTIONS

- Do not spray directly into circulatory pathways. Not for intravascular application.
- Application of the liquids must be performed under visual control.
- Single use only. Do not attempt to re-sterilize or reuse device.
- Device may degrade if reprocessed. No effective cleaning process has been developed to allow reuse of device. Reuse of device may lead to injury or illness of a patient.

For safe and proper use of the GASLESS TIP device, please refer to the full Instructions for Use.

For questions or ordering information, please contact your
Baxter representative or call 888-229-0001.

www.advancedsurgery.baxter.com

Advancing the art of healing

USA
Baxter Healthcare Corporation
1 Baxter Parkway
Deerfield, IL 60015

1. Gasless Tip Applicator Instructions for Use

Baxter, Advancing the art of healing and Tisseel are trademarks of Baxter International Inc.
©2020 Baxter Healthcare Corporation. All rights reserved.

USMP/90/20-0031 5/2020

Baxter